

Microscopes



GrayMed
 PH: 61 3 86091732
 Email: sales@graymed.com.au
 Web: www.graymed.com.au

Click on microscope picture for more information



| | | | | |
|---------------------|-----------------------|--|--|--|
| BEST SELLERS | Model | Mi5 | i4 Labscope | Revelation III |
| | Model # | M5M-B05A-IPL3 | I4M-BN4A-IPL3 | R3M-BN4A-DAL3 |
| | Head | Binoc | Binoc | Binoc |
| | Objectives | 4-10-20-40-100 Infinity Plan | 4-10-40-100 Infinity Plan | 4-10-40-100 DIN Achromat |
| | Nosepiece | Reverse quintuple | Reverse quadruple | Quadruple forward-facing |
| | Stage | Mechanical | Mechanical | Mechanical |
| | Target Markets | Pathology, Hematology, Oncology, Gynecology, Research, Microbiology, IHC | Cytology, Hematology, Oncology, Gynecology, Research, Microbiology | Routine Lab, POL, Tech Schools, Universities |

Click on microscope picture for more information



| | | | | |
|--------------------|-----------------------|---|---|----------------------------------|
| SPECIALIZED | Model | Mi5 Polarizer | Lumin™ i4 | i4 Semen Evaluation |
| | Model # | M5M-BPOL-IPH3 | I4S-EPB4-IPL3 | i4S-SEB4-iPL3 |
| | Head | Binoc | Binoc | Binoc |
| | Objectives | 4-10-20-40 Infinity Plan | 4-10-40-100 Infinity Plan | 4-10-40-100 DIN Phase Plan |
| | Nosepiece | Reverse quintuple | Reverse quadruple | Reverse quadruple |
| | Stage | Round, Rotates 360, Mechanical Optional | Mechanical | Mechanical |
| | Target Markets | Rheumatology, Urology | Tuberculosis, Malaria, STD, Pathology, Immunology, Research IHC | Fertility Clinics, IVF, Research |



P 770.270.1394
 F 770.270.2389

865 Marathon Parkway
 Lawrenceville GA 30046

LWScientific.com