

# Centrifuges



GrayMed  
 PH: 61 3 86091732  
 Email: sales@graymed.com.au  
 Web: www.graymed.com.au



Click on centrifuge picture for more information

Model	Universal Fixed	*E8 Digital	*E8 Variable	*E8 Fixed
Model #	UNC-08AD-15T3	E8C-U8AD-15T3	E8C-U8AV-1503	E8C-U8AF-1503
Type	Test Tube	Test Tube	Test Tube	Test Tube
Features	Programmable	Programmable	Variable Speed	Fixed Speed
Capacity	8 place X 15ml	8 place X 15ml	8 place X 15ml	8 place X 15ml
Tube Size Min/Max	10.25 X 65mm / 17 X 130mm	10.25 X 65mm / 17 X 130mm	10.25 X 65mm / 17 X 130mm	10.25 X 65mm / 17 X 130mm
Max g	1350	1400	1400	1500
Max Speed (rpm)	3300	+/-3500	+/-3500	+/-3500
Applications	Blood, Urine, Semen, Fecal	Blood, Urine, Semen, Fecal	Blood, Urine, Semen, Fecal	Blood Only

\*12V Portafuge option available (includes carseat)

Click on centrifuge picture for more information



Model	Universal Swing-Out	C5	C5 - 24 Place
Model #	UNC-06SD-15T3	C5C-08SU-15T3	C5C-BKSU-15T3
Type	Test Tube	Test Tube Swing Out	Test Tube Swing Out
Features	Programmable	5000 RPM Programmable	5000 RPM Programmable
Capacity	6 place X 15ml	8 place X 15ml *4 place X 50ml	24 X 2.5-10ml *4 place X 50ml tube inserts
Tube Size Min/Max	10.25 X 65mm / 17 X 130mm	10.25 X 65mm / 16 X 130mm	10.25 X 65mm / 17 X 130mm
Max g	1350	4000	4000
Max Speed (rpm)	3300	5000	5000
Applications	Blood, Urine, Semen, Fecal	Blood (STAT), Urine, Semen, Fecal, Hematocrits, PRP/PPP, Stem-Cell, Irap	Blood (STAT), Urine, Semen, Fecal, Hematocrits, PRP/PPP, Stem-Cell

\*Optional Features



P 770.270.1394  
 F 770.270.2389

865 Marathon Parkway  
 Lawrenceville GA 30046

LWScientific.com

# Centrifuges



GrayMed  
 PH: 61 3 86091732  
 Email: sales@graymed.com.au  
 Web: www.graymed.com.au



Click on centrifuge picture for more information

Model	V24 Combo	V24 Combo Tach	M24 Hematocrit
Model #	CMP-BX7N-7701 (*choose rotor)	CMP-BX7N-77T1 (*choose rotor)	HEC-24HF-7501
Type	Combination	Combination	Microhematocrit
Features	Test Tube & Microhematocrit	Test Tube & Microhematocrit	Heavy Duty 24 Place
Capacity	24 X 15ml or 24 X 75mm	24 X 15ml or 24 X 75mm	24 X 75mm
Tube Size Min/Max	10.25 X 65mm /17X130mm or 75mm capillary	10.25X65mm/17X130mm or 75mm capillary	75mm capillary (only)
Max g	2450 - 14900	2450 - 14900	14900
Max Speed (rpm)	4000 - 12000	4000 - 12000	4000 - 12000
Applications	Blood, Urine, Semen, Fecal, Hematocrits, DNA Separation	Blood, Urine, Semen, Fecal, Hematocrits, DNA Separation	Hematocrits Only

**\*Optional Rotors for Combo V24**



24-Place Rotor (CMP-RT24-15TT)



12-Place Rotor (CMP-RT12-15TT)



24-Place Microtube Rotor (CMP-RT24-02MT)



24-Place Hematocrit Rotor (CMP-RT24-75HE)

Click on centrifuge picture for more information



Model	Zipocrit	Digital ZipSpin
Model #	ZPC-04HF-7503	ZPC-08AD-0203
Type	Microhematocrit	Microcentrifuge
Features	Portable angled hematocrit	8-Place Microtube or 16-Place PCR
Capacity	4 X 75mm	1.5 -2.0ml (8-place) 2x8x0.2ml (16-place)
Tube Size Min/Max	75mm capillary (only)	0.5ml / 2.0ml
Max g	4400	2000
Max Speed (rpm)	10000	1000-6000
Applications	Hematocrits Only	Blood, Urine, DNA Separation, HPLC microfiltration



P 770.270.1394  
 F 770.270.2389

865 Marathon Parkway  
 Lawrenceville GA 30046

LWScientific.com