

CENTRIFUGE COMPARISON CHART

SWING-OUT

ANGLED



	Combo XL	Universal	E8 Digital	E8 Variable	E8 Fixed	Zip-IQ TT	Universal
MODEL	Combo XL	Universal	E8 Digital	E8 Variable	E8 Fixed	Zip-IQ TT	Universal
MODEL NUMBER	CMC-16SD-50T3	UNC-06SD-15T3	E8C-U8AD-15T3	E8C-U8AV-1503	E8C-U8AF-1503	ZIC-06AD-15T3	UNC-08AD-15T3
TYPE	Test Tube	Test Tube	Test Tube	Test Tube	Test Tube	Test Tube	Test Tube
FEATURES	Programmable	Programmable	Programmable	Variable Speed	Fixed Speed	Programmable	Programmable
CAPACITY	12 place X 15ml & 4 place 50ml	6 place X 15ml	8 place X 15ml	8 place X 15ml	8 place X 15ml	6 place X 15ml	8 place X 15ml
MAX G-FORCE	2683	1800	1500	1500	1500	2500	1450
MAX SPEED	4000	3400	3500	3500	3500	5000	3400
TUBE SIZE (MIN/MAX)	10.25 X 65mm / 17 X 130mm / 29.5 X 118mm	10.25 X 65mm / 17 X 130mm	10.25 X 65mm / 17 X 130mm	10.25 X 65mm / 17 X 130mm	10.25 X 65mm / 17 X 130mm	13 x 75mm / 16 x 100mm	10.25 X 65mm / 17 X 130mm
APPLICATIONS	Blood, Urine, Semen, Fecal	Blood, Urine, Semen, Fecal	Blood, Urine, Semen, Fecal	Blood, Urine, Semen, Fecal	Blood Only	Blood (STAT), Urine, Semen, Fecal	Blood, Urine, Semen, Fecal

CENTRIFUGE COMPARISON CHART

HEMATOCRIT

COMBINATION

PORTABLE

HEATED



	Zip-IQ PCV	V24 Combo w/ Tach	V24 Combo	ZipCombo	Zip-IQ MT / PCR	Portafuge	Universal
MODEL	Zip-IQ PCV	V24 Combo w/ Tach	V24 Combo	ZipCombo	Zip-IQ MT / PCR	Portafuge	Universal
MODEL NUMBER	ZIC-24HD-75T3	CMP-BX7N-77T1	CMP-BX7N-7701	ZCC-1206-77T3	ZiC-24AD-02T3	E8C-U8AF-150P*	PTC-06SH-15T3
TYPE	Hematocrit	Combination*	Combination*	Combination	Combination	Test Tube	Test Tube
FEATURES	Programmable	Test Tube & Hematocrit	Test Tube & Hematocrit	Microtube & Hematocrit	Microtube & PCR Strips	Portable	Heated Bowl Chamber: ambient to 71°C (160°F)
CAPACITY	24 place X 75mm	24 X 15ml or 24 X 75mm	24 X 15ml or 24 X 75mm	12 place X 40mm / 6 place 0.2-2.0ml	12 X 0.5-0.9ml; 12 X 1.0-2.0ml / 4 x 8-place PCR Strips	8 place X 15ml	6 place X 15ml
MAX G-FORCE	13680	14900	14900	6900	4850	1534	759
MAX SPEED	12000	12000	12000	12000	8500	3500	2200
TUBE SIZE (MIN/MAX)	75mm capillary (only)	10.25 X 65mm / 17 X 130mm or 75mm capillary	10.25 X 65mm / 17 X 130mm or 75mm capillary	40 X 0.95mm / 0.5ml - 2.0ml	0.5ml / 2.0ml	10.25 X 65mm / 17 X 130mm	10.25 X 65mm / 17 X 130mm
APPLICATIONS	Hematocrits Only	Blood, Urine, Semen, Fecal, Hematocrits, DNA Separation	Blood, Urine, Semen, Fecal, Hematocrits, DNA Separation	Blood, Urine, DNA Separation, HPLC Microfiltration	Blood, Urine, DNA Separation, HPLC Microfiltration	Blood Only	Oil Testing

*4 rotor choices available

*4 rotor choices available

*Variable and Digital models available