

# Zipocrit Centrifuge



GrayMed  
PH: 61 3 86091732  
Email: sales@graymed.com.au  
Web: www.graymed.com.au

## SPECIFICATIONS

- Portable and compact microhematocrit centrifuge
- Perfect for medical missionaries or veterinary field work
- Spins four 75mm capillary tubes in angled rotor
- Mylar-coated or plastic capillary tubes recommended - vial of 200 tubes included; Critoseal clay available
- Speed: 10,000 rpm
- RCF: 4,400g
- Fifteen minute timer - separations in five minutes
- Designed with a small footprint to conserve space
- Safety switch lid
- Strobe through transparent lid for exact speed
- LED power indicator light
- DC 12V power (AC adapter power cord included)
- Optional 12V DC car battery cables and 220V AC adapter power cord available



## Zipocrit Centrifuge

Model #	Description
ZPC-04HF-7503	Zipocrit centrifuge

## DIMENSIONS AND WEIGHT

Height: 4.5" (114.3mm)  
Length: 8" (203mm)  
Width: 6" (152mm)  
Weight: 3.3lbs. (1.5kg)



GrayMed  
PH: 61 3 86091732  
Email: sales@graymed.com.au  
Web: www.graymed.com.au



P 770.270.1394  
F 770.270.2389

865 Marathon Parkway  
Lawrenceville GA 30046

LWScientific.com